



NTF2 Antibody

Product Code	CSB-PA016214XA01CZF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P926
Immunogen	Recombinant Candida albicans (Yeast) NTF2 protein
Raised In	Rabbit
Species Reactivity	Candida albicans (Yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Candida albicans (Yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	NTF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.