



NUDT3 Antibody

Product Code	CSB-PA530966XA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95989
Immunogen	Recombinant Homo sapiens (Human) NUDT3 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens (Human)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (14-16 weeks)
Target Names	NUDT3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.