



PHO85 Antibody

Product Code	CSB-PA875759XA01CZF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HGY5
Immunogen	Recombinant Candida albicans (Yeast) PHO85 protein
Raised In	Rabbit
Species Reactivity	Candida albicans (Yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Candida albicans (Yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	PHO85
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.