



PRM5 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA522241XA01SVO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | E7KPQ3 |
| Immunogen | Recombinant Saccharomyces cerevisiae (strain Lalvin QA23) (Baker's yeast) PRM5 protein |
| Raised In | Rabbit |
| Species Reactivity | Saccharomyces cerevisiae (strain Lalvin QA23) (Baker's yeast) |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Saccharomyces cerevisiae (strain Lalvin QA23) (Baker's yeast) |
| Lead Time | Made-to-order (14-16 weeks) |
| Target Names | PRM5 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |