



# Recombinant *Macaca fascicularis* Tubulin alpha-4A chain (TUBA4A)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-MP025313MOV  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.  |
| <b>Uniprot No.</b>       | P68367   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | <i>Macaca fascicularis</i> (Crab-eating macaque) ( <i>Cynomolgus</i> monkey)   |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MRECISVHVG QAGVQMGNAC WELYCLEHGI QPDGQMPSDK<br>TIGGGDSSFT TFFCETGAGK HVPRAVFVDL EPTVIDEIRN GPYRQLFHPE<br>QLITGKEDAA NNYARGHYTI GKEIIDPVLDRIRKLSLSDQCT GLQGFLVFHS<br>FGGGTGSFT SLLMERLSVD YGKKSLEFS IYPAPQVSTA VVEPYNSILT<br>THTTLEHSDC AFMVDNEAIY DICRRNLIDIE RPTYTNLNLRL ISQIVSSITA<br>SLRFDGALNV DLTEFQTNLV PYPRIHFPLA TYAPVISA EK AYHEQLSVAE<br>ITNACFEPAN QMVKCDPRHG KYMACCLLYR GDVVPKDVNA AIAAIKTKRS<br>IQFVDWCPTG FKVGINYQPP TVVPGDLAK VQRAVCMLSN TTAIAEAWAR<br>LDHKFDL MYA KRAVHWYVG EGMEEGEFSE AREDMAALEK DYEEVGIDSY<br>EDEDEGEE |
| <b>Source</b>            | Mammalian cell   |
| <b>Target Names</b>      | TUBA4A   |
| <b>Protein Names</b>     | Recommended name: Tubulin alpha-4A chain Alternative name(s): Alpha-tubulin 1 Testis-specific alpha-tubulin Tubulin alpha-1 chain  |
| <b>Expression Region</b> | 1-448  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |