



URA9 Antibody

Product Code	CSB-PA764888XA01KAN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6SZS6
Immunogen	Recombinant Kluyveromyces marxianus (Yeast) (Candida kefyr) URA9 protein
Raised In	Rabbit
Species Reactivity	Kluyveromyces marxianus (Yeast) (Candida kefyr)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Kluyveromyces marxianus (Yeast) (Candida kefyr)
Lead Time	Made-to-order (14-16 weeks)
Target Names	URA9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.