



RIB7 Antibody

Product Code	CSB-PA865225XA01KAN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P4B8
Immunogen	Recombinant Kluyveromyces marxianus (Yeast) (Candida kefyr) RIB7 protein
Raised In	Rabbit
Species Reactivity	Kluyveromyces marxianus (Yeast) (Candida kefyr)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Kluyveromyces marxianus (Yeast) (Candida kefyr)
Lead Time	Made-to-order (14-16 weeks)
Target Names	RIB7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.