



RKM5 Antibody

Product Code	CSB-PA517213XA01SVQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E7LXI9
Immunogen	Recombinant Saccharomyces cerevisiae (strain VIN 13) (Baker's yeast) RKM5 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae (strain VIN 13) (Baker's yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces cerevisiae (strain VIN 13) (Baker's yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	RKM5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.