



RPL10A Antibody

Product Code	CSB-PA886801XA01DSJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9SW75
Immunogen	Recombinant Chlamydomonas reinhardtii (Chlamydomonas smithii) RPL10A protein
Raised In	Rabbit
Species Reactivity	Chlamydomonas reinhardtii (Chlamydomonas smithii)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Chlamydomonas reinhardtii (Chlamydomonas smithii)
Lead Time	Made-to-order (14-16 weeks)
Target Names	RPL10A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.