



RRG8 Antibody

Product Code	CSB-PA477655XA01SVH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B5VTQ6
Immunogen	Recombinant Saccharomyces cerevisiae (strain AWRI1631) (Baker's yeast) RRG8 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae (strain AWRI1631) (Baker's yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces cerevisiae (strain AWRI1631) (Baker's yeast)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	RRG8