



POLR3G Antibody

Product Code	CSB-PA018349XA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15318
Immunogen	Recombinant Homo sapiens (Human) POLR3G protein
Raised In	Rabbit
Species Reactivity	Homo sapiens (Human)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (14-16 weeks)
Target Names	POLR3G
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.