



RPS1B Antibody

Product Code	CSB-PA512057XA01SVL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C7GVZ9
Immunogen	Recombinant Saccharomyces cerevisiae (strain JAY291) (Baker's yeast) RPS1B protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae (strain JAY291) (Baker's yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces cerevisiae (strain JAY291) (Baker's yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	RPS1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.