



RTC4 Antibody

Product Code	CSB-PA454237XA01ZBK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B2G493
Immunogen	Recombinant Zygosaccharomyces rouxii (Candida mogii) RTC4 protein
Raised In	Rabbit
Species Reactivity	Zygosaccharomyces rouxii (Candida mogii)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Zygosaccharomyces rouxii (Candida mogii)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	RTC4