



# SAP190 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA514968XA01SVN  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | C8ZCJ2   |
| <b>Immunogen</b>           | Recombinant Saccharomyces cerevisiae (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast) SAP190 protein |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Saccharomyces cerevisiae (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast)                            |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)   |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                             |
| <b>Purification Method</b> | Antigen Affinity Purified  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Product Type</b>        | Polyclonal Antibody  |
| <b>Immunogen Species</b>   | Saccharomyces cerevisiae (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast)                            |
| <b>Lead Time</b>           | Made-to-order (14-16 weeks)  |
| <b>Target Names</b>        | SAP190   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                                  |