



# YIL171W-A Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA208134XA01SVG   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                       |
| <b>Uniprot No.</b>         | A0A023PXM9  |
| <b>Immunogen</b>           | Recombinant Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast) YIL171W-A protein |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)                               |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)  |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Non-conjugated  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                    |
| <b>Purification Method</b> | Antigen Affinity Purified   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Product Type</b>        | Polyclonal Antibody   |
| <b>Immunogen Species</b>   | Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)                               |
| <b>Lead Time</b>           | Made-to-order (14-16 weeks)   |
| <b>Target Names</b>        | YIL171W-A   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                         |