



EC1118_1M3_2520g Antibody

| | |
|----------------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA513912XA01SVN |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | C8ZET7 |
| Immunogen | Recombinant Saccharomyces cerevisiae (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast) EC1118_1M3_2520g protein |
| Raised In | Rabbit |
| Species Reactivity | Saccharomyces cerevisiae (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast) |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Saccharomyces cerevisiae (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast) |
| Lead Time | Made-to-order (14-16 weeks) |
| Gene Names | EC1118_1M3_2520g |