



# YER023C-A Antibody

|                     |   |
|---------------------|---|
| Product Code        | CSB-PA848437XA01SVG   |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                       |
| Uniprot No.         | Q8TGP2  |
| Immunogen           | Recombinant Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast) YER023C-A protein |
| Raised In           | Rabbit  |
| Species Reactivity  | Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)                               |
| Tested Applications | ELISA, WB (ensure identification of antigen)  |
| Form                | Liquid  |
| Conjugate           | Non-conjugated  |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                    |
| Purification Method | Antigen Affinity Purified   |
| Isotype             | IgG   |
| Clonality           | Polyclonal  |
| Product Type        | Polyclonal Antibody   |
| Immunogen Species   | Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)                               |
| Lead Time           | Made-to-order (14-16 weeks)   |
| Gene Names          | YER023C-A   |