



EC1118_1L10_0100g Antibody

Product Code	CSB-PA513181XA01SVN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C8ZCS2
Immunogen	Recombinant Saccharomyces cerevisiae (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast) EC1118_1L10_0100g protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces cerevisiae (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	EC1118_1L10_0100g
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.