



SKN1 Antibody

Product Code	CSB-PA164227XA01CZF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P87024
Immunogen	Recombinant Candida albicans (Yeast) SKN1 protein
Raised In	Rabbit
Species Reactivity	Candida albicans (Yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Candida albicans (Yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	SKN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.