



# SLX1 Antibody

|                            |                                                                                                                                     |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA751803XA01DOT                                                                                                                 |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                       |
| <b>Uniprot No.</b>         | Q75EA5                                                                                                                              |
| <b>Immunogen</b>           | Recombinant Ashbya gossypii (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (Eremothecium gossypii) SLX1 protein |
| <b>Raised In</b>           | Rabbit                                                                                                                              |
| <b>Species Reactivity</b>  | Ashbya gossypii (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (Eremothecium gossypii)                          |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)                                                                                        |
| <b>Form</b>                | Liquid                                                                                                                              |
| <b>Conjugate</b>           | Non-conjugated                                                                                                                      |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                                                    |
| <b>Purification Method</b> | Antigen Affinity Purified                                                                                                           |
| <b>Isotype</b>             | IgG                                                                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                                                                          |
| <b>Product Type</b>        | Polyclonal Antibody                                                                                                                 |
| <b>Immunogen Species</b>   | Ashbya gossypii (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (Eremothecium gossypii)                          |
| <b>Lead Time</b>           | Made-to-order (14-16 weeks)                                                                                                         |
| <b>Target Names</b>        | SLX1                                                                                                                                |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                                                         |