



SHO1 Antibody

Product Code	CSB-PA879146XA01EPX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P864
Immunogen	Recombinant Cyberlindnera jadinii (Torula yeast) (Pichia jadinii) SHO1 protein
Raised In	Rabbit
Species Reactivity	Cyberlindnera jadinii (Torula yeast) (Pichia jadinii)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Cyberlindnera jadinii (Torula yeast) (Pichia jadinii)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	SHO1