



SNF1 Antibody

Product Code	CSB-PA345995XA01CZF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52497
Immunogen	Recombinant Candida albicans (Yeast) SNF1 protein
Raised In	Rabbit
Species Reactivity	Candida albicans (Yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Candida albicans (Yeast)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	SNF1