



TPP1 Antibody

Product Code	CSB-PA864050XA01GHS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HET9
Immunogen	Recombinant Saccharomyces mikatae (Yeast) TPP1 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces mikatae (Yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces mikatae (Yeast)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	TPP1