



TDA7 Antibody

Product Code	CSB-PA518968XA01SVR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E7QJS8
Immunogen	Recombinant Saccharomyces cerevisiae (strain Zymaflore VL3) (Baker's yeast) TDA7 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae (strain Zymaflore VL3) (Baker's yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces cerevisiae (strain Zymaflore VL3) (Baker's yeast)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	TDA7