



ugt52 Antibody

Product Code	CSB-PA706744XA01DKK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q54IL5
Immunogen	Recombinant Dictyostelium discoideum (Slime mold) ugt52 protein
Raised In	Rabbit
Species Reactivity	Dictyostelium discoideum (Slime mold)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Dictyostelium discoideum (Slime mold)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ugt52
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.