



URE2 Antibody

Product Code	CSB-PA759178XA01FZP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7LLZ8
Immunogen	Recombinant Saccharomyces paradoxus (Yeast) (Saccharomyces douglasii) URE2 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces paradoxus (Yeast) (Saccharomyces douglasii)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces paradoxus (Yeast) (Saccharomyces douglasii)
Lead Time	Made-to-order (14-16 weeks)
Target Names	URE2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.