



E2 Antibody

Product Code	CSB-PA767977XA01HOM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q84294
Immunogen	Recombinant Human papillomavirus type 6a E2 protein
Raised In	Rabbit
Species Reactivity	Human papillomavirus type 6a
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Human papillomavirus type 6a
Lead Time	Made-to-order (14-16 weeks)
Target Names	E2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.