



xynA Antibody

Product Code	CSB-PA299151XA01AQY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12562
Immunogen	Recombinant Aureobasidium pullulans (Black yeast) (Pullularia pullulans) xynA protein
Raised In	Rabbit
Species Reactivity	Aureobasidium pullulans (Black yeast) (Pullularia pullulans)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Aureobasidium pullulans (Black yeast) (Pullularia pullulans)
Lead Time	Made-to-order (14-16 weeks)
Target Names	xynA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.