



Xyn2 Antibody

Product Code	CSB-PA549916XA01TKN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B5A7N4
Immunogen	Recombinant Trichoderma harzianum (Hypocrea lixii) Xyn2 protein
Raised In	Rabbit
Species Reactivity	Trichoderma harzianum (Hypocrea lixii)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Trichoderma harzianum (Hypocrea lixii)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Xyn2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.