



xyn2 Antibody

Product Code	CSB-PA281460XA01GHW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	G0RUP7
Immunogen	Recombinant Hypocrea jecorina (strain QM6a) (Trichoderma reesei) xyn2 protein
Raised In	Rabbit
Species Reactivity	Hypocrea jecorina (strain QM6a) (Trichoderma reesei)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Hypocrea jecorina (strain QM6a) (Trichoderma reesei)
Lead Time	Made-to-order (14-16 weeks)
Target Names	xyn2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.