



L1 Antibody

Product Code	CSB-PA300898XA01HOM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69898
Immunogen	Recombinant Human papillomavirus type 6a L1 protein
Raised In	Rabbit
Species Reactivity	Human papillomavirus type 6a
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Human papillomavirus type 6a
Lead Time	Made-to-order (14-16 weeks)
Target Names	L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.