



ZRG8 Antibody

Product Code	CSB-PA408133XA01STA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6ZQY2
Immunogen	Recombinant Saccharomyces cerevisiae (strain YJM789) (Baker's yeast) ZRG8 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae (strain YJM789) (Baker's yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces cerevisiae (strain YJM789) (Baker's yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ZRG8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.