



TY1A-PR3 Antibody

Product Code	CSB-PA750902XA01SVG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6Q5H1
Immunogen	Recombinant Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast) TY1A-PR3 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	TY1A-PR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.