



# YIR030W-A Antibody

|                     |                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------|
| Product Code        | CSB-PA166521XA01SVG                                                                                 |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                       |
| Uniprot No.         | A0A023PZK4                                                                                          |
| Immunogen           | Recombinant Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast) YIR030W-A protein |
| Raised In           | Rabbit                                                                                              |
| Species Reactivity  | Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)                               |
| Tested Applications | ELISA, WB (ensure identification of antigen)                                                        |
| Form                | Liquid                                                                                              |
| Conjugate           | Non-conjugated                                                                                      |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                    |
| Purification Method | Antigen Affinity Purified                                                                           |
| Isotype             | IgG                                                                                                 |
| Clonality           | Polyclonal                                                                                          |
| Product Type        | Polyclonal Antibody                                                                                 |
| Immunogen Species   | Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)                               |
| Lead Time           | Made-to-order (14-16 weeks)                                                                         |
| Target Names        | YIR030W-A                                                                                           |
| Usage               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                         |