



AOX1 Antibody

Product Code	CSB-PA322107XA01GHX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22185
Immunogen	Recombinant Sauromatum venosum (Voodoo lily) (Typhonium venosum) AOX1 protein
Raised In	Rabbit
Species Reactivity	Sauromatum venosum (Voodoo lily) (Typhonium venosum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Sauromatum venosum (Voodoo lily) (Typhonium venosum)
Lead Time	Made-to-order (14-16 weeks)
Target Names	AOX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.