



# SPCC23B6.02c Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA891873XA01SXV  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | Q9UUA0   |
| <b>Immunogen</b>           | Recombinant Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) SPCC23B6.02c protein |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)                                  |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)   |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                     |
| <b>Purification Method</b> | Antigen Affinity Purified  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Product Type</b>        | Polyclonal Antibody  |
| <b>Immunogen Species</b>   | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)                                  |
| <b>Lead Time</b>           | Made-to-order (14-16 weeks)  |
| <b>Target Names</b>        | SPCC23B6.02c   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                          |