



matMi Antibody

Product Code	CSB-PA744302XA01GHY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6WRX7
Immunogen	Recombinant Schizosaccharomyces kambucha (Fission yeast) matMi protein
Raised In	Rabbit
Species Reactivity	Schizosaccharomyces kambucha (Fission yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Schizosaccharomyces kambucha (Fission yeast)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	matMi