



SPAC2E1P5.01c Antibody

Product Code	CSB-PA879117XA01SXV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P7C3
Immunogen	Recombinant Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) SPAC2E1P5.01c protein
Raised In	Rabbit
Species Reactivity	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	SPAC2E1P5.01c
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.