



# nce101 Antibody

|                            |                                                                                                |
|----------------------------|------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA511303XA01SXV                                                                            |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                  |
| <b>Uniprot No.</b>         | C6Y4B6                                                                                         |
| <b>Immunogen</b>           | Recombinant Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) nce101 protein |
| <b>Raised In</b>           | Rabbit                                                                                         |
| <b>Species Reactivity</b>  | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)                            |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)                                                   |
| <b>Form</b>                | Liquid                                                                                         |
| <b>Conjugate</b>           | Non-conjugated                                                                                 |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4               |
| <b>Purification Method</b> | Antigen Affinity Purified                                                                      |
| <b>Isotype</b>             | IgG                                                                                            |
| <b>Clonality</b>           | Polyclonal                                                                                     |
| <b>Product Type</b>        | Polyclonal Antibody                                                                            |
| <b>Immunogen Species</b>   | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)                            |
| <b>Lead Time</b>           | Made-to-order (14-16 weeks)                                                                    |
| <b>Gene Names</b>          | nce101                                                                                         |