



RAD17 Antibody

Product Code	CSB-PA894555XA01DSU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9XT62
Immunogen	Recombinant Chlorocebus aethiops (Green monkey) (Cercopithecus aethiops) RAD17 protein
Raised In	Rabbit
Species Reactivity	Chlorocebus aethiops (Green monkey) (Cercopithecus aethiops)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Chlorocebus aethiops (Green monkey) (Cercopithecus aethiops)
Lead Time	Made-to-order (14-16 weeks)
Target Names	RAD17
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.