



RAN1 Antibody

Product Code	CSB-PA339981XA01SZY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38546
Immunogen	Recombinant Solanum lycopersicum (Tomato) (Lycopersicon esculentum) RAN1 protein
Raised In	Rabbit
Species Reactivity	Solanum lycopersicum (Tomato) (Lycopersicon esculentum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Solanum lycopersicum (Tomato) (Lycopersicon esculentum)
Lead Time	Made-to-order (14-16 weeks)
Target Names	RAN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.