



rad9 Antibody

Product Code	CSB-PA342919XA01GHZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48013
Immunogen	Recombinant Schizosaccharomyces octosporus (Fission yeast) (Octosporomyces octosporus) rad9 protein
Raised In	Rabbit
Species Reactivity	Schizosaccharomyces octosporus (Fission yeast) (Octosporomyces octosporus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Schizosaccharomyces octosporus (Fission yeast) (Octosporomyces octosporus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	rad9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.