



# SPAC19D5.10c Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA737606XA01SXV  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | Q65ZA5   |
| <b>Immunogen</b>           | Recombinant Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) SPAC19D5.10c protein |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)                                  |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)   |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                     |
| <b>Purification Method</b> | Antigen Affinity Purified  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Product Type</b>        | Polyclonal Antibody  |
| <b>Immunogen Species</b>   | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)                                  |
| <b>Lead Time</b>           | Made-to-order (14-16 weeks)  |
| <b>Gene Names</b>          | SPAC19D5.10c   |