



SPAC6F6.19 Antibody

Product Code	CSB-PA514900XA01SXV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C6Y4A5
Immunogen	Recombinant Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) SPAC6F6.19 protein
Raised In	Rabbit
Species Reactivity	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	SPAC6F6.19
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.