



# Recombinant Danio rerio Ribonucleoside-diphosphate reductase subunit M2 (rrm2)

<b>Product Code</b>	CSB-BP020519DIL
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P79733
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSSTRSPLKT KNENTISTKM NNMSFVDKEN TPPSLSSSTRI LASKTARKIF DESEGQSKAK KGAVEEEPLL KENPHRFVIF PIQYHDIWQM YKKAESFWT AEEVDLSKDL QHWDSLKDEE RYFISHVLAF FAASDGVNE NLVERFTQEV QVTEARCFYG FQIAMENIHS EMYSLIDTY IKDSKEREFL FNAIETMPCV KKKADWALNW IGDKNARYGE RVVAFAAVEG IFFSGSFASI FWLKKRGLMP GLTFSNELIS RDEGLHCDFA CLMFKHLINK PSEETVKKII MNAVEIEQEF LTDALPVKLI GMNCDLMKQY IEFVADRLLL ELGFDKVYRV ENPFDPMENI SLEGKTNFFE KRVGEYQRMG VMSGTTDNTF TLDADF
<b>Source</b>	Baculovirus
<b>Target Names</b>	rrm2
<b>Protein Names</b>	Recommended name: Ribonucleoside-diphosphate reductase subunit M2 EC=1.17.4.1 Alternative name(s): Ribonucleotide reductase protein R2 class I Ribonucleotide reductase small chain Ribonucleotide reductase small subunit
<b>Expression Region</b>	1-386
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.