



# SPAC27E2.11c Antibody

|                            |                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA892223XA01SXV                                                                                  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                        |
| <b>Uniprot No.</b>         | Q9UTA0                                                                                               |
| <b>Immunogen</b>           | Recombinant Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) SPAC27E2.11c protein |
| <b>Raised In</b>           | Rabbit                                                                                               |
| <b>Species Reactivity</b>  | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)                                  |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)                                                         |
| <b>Form</b>                | Liquid                                                                                               |
| <b>Conjugate</b>           | Non-conjugated                                                                                       |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                     |
| <b>Purification Method</b> | Antigen Affinity Purified                                                                            |
| <b>Isotype</b>             | IgG                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                           |
| <b>Product Type</b>        | Polyclonal Antibody                                                                                  |
| <b>Immunogen Species</b>   | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)                                  |
| <b>Lead Time</b>           | Made-to-order (14-16 weeks)                                                                          |
| <b>Target Names</b>        | SPAC27E2.11c                                                                                         |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                          |