



Recombinant 6-phosphofructokinase (pfkA)

Product Code	CSB-MP017822LMM
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P80019
Product Type	Recombinant Protein
Immunogen Species	Lactobacillus delbrueckii subsp. bulgaricus
Purity	≥85% (SDS-PAGE)
Sequence	MKRIGILTSG GDAPGMNAAV RAVTRVAIAN GLEVFGIRYG FAGLVAGDIF PLESEDVAHL INVSQTFLYS ARYPEFAEEE GQLAGIEQLK KHGIDAVVVI GGDGSYHGAL QLTRHGFNSI GLPGTIDNDI PYTDATIGYD TACMTAMDAI DKIRDTASSH HRVFIVNVMG RNCGDIA MRV GVACGADAIV IPERPYDVEE IANRLKQAQE SGKDHGLVVV AEGVMTADQF MAELKKYGDF DVRANVLGHM QRGGTPTVSD RVLASKLGSE AVHLLLEGKG GLAVGIENGK VTSHDILDLF DESHRGDYDL LKLNADLSR
Source	Mammalian cell
Target Names	pfkA
Protein Names	Recommended name: 6-phosphofructokinase Short name= Phosphofructokinase EC= 2.7.1.11 Alternative name(s): Phosphohexokinase
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.