



# SPCC569.06 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA897341XA01SXV  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                      |
| <b>Uniprot No.</b>         | Q9Y7S5   |
| <b>Immunogen</b>           | Recombinant Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) SPCC569.06 protein |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)                                |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)   |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                   |
| <b>Purification Method</b> | Antigen Affinity Purified  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Product Type</b>        | Polyclonal Antibody  |
| <b>Immunogen Species</b>   | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)                                |
| <b>Lead Time</b>           | Made-to-order (14-16 weeks)  |
| <b>Gene Names</b>          | SPCC569.06   |