



# SPBC83.12 Antibody

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| Product Code        | CSB-PA529727XA01SXV                                                                               |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                     |
| Uniprot No.         | O94696                                                                                            |
| Immunogen           | Recombinant Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) SPBC83.12 protein |
| Raised In           | Rabbit                                                                                            |
| Species Reactivity  | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)                               |
| Tested Applications | ELISA, WB (ensure identification of antigen)                                                      |
| Form                | Liquid                                                                                            |
| Conjugate           | Non-conjugated                                                                                    |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                  |
| Purification Method | Antigen Affinity Purified                                                                         |
| Isotype             | IgG                                                                                               |
| Clonality           | Polyclonal                                                                                        |
| Product Type        | Polyclonal Antibody                                                                               |
| Immunogen Species   | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)                               |
| Lead Time           | Made-to-order (14-16 weeks)                                                                       |
| Gene Names          | SPBC83.12                                                                                         |