



SPBC56F2.05c Antibody

| | |
|----------------------------|------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA526431XA01SXV |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O60056 |
| Immunogen | Recombinant Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) SPBC56F2.05c protein |
| Raised In | Rabbit |
| Species Reactivity | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) |
| Lead Time | Made-to-order (14-16 weeks) |
| Target Names | SPBC56F2.05c |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |